
## Instrumententische ST Serie/ Instrument tables ST series

## Gebrauchsanweisung/ direction for use Version 2.2

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Before using the devices, please make yourself familiar with the instructions and safety guidelines. Read the instructions carefully and follow the notes on safe handling of the device.

We reserve the right to change design and delivery due to technical advancements.

**Please note**

**Table of contents**

Instructions .....	35
Medical directive.....	35
Transport regulations.....	36
Warning and information signs .....	37
General construction .....	38
Operation.....	40
Intended use.....	41
Instructions for assembly and usage .....	42
Change of main fuses.....	45
Variations of system tables.....	47
Patient surrounding .....	52
Instrument table with isolation transformer .....	53
Disposal.....	55
Behavior in case of failure .....	57
Environmental conditions .....	58
Electromagnetic Emissions .....	59
Order data - spare parts .....	64

## Instructions

- Instrument tables ST series were constructed regarding the Standard DIN EN 60601-1. Manufacturing, testing, installation, maintenance and repair occur in compliance with German and international requirements.
- Do not use as a means of transport
- All electrical components connected to line voltage comply with VDE.
- The instrument tables of the ST series are in accordance with the EU Directive 2017/745 (MDR) on Medical Devices for active devices of class I - non-invasive.
- Accident prevention measures due to legal requirements shall be followed by the user.
- The instrument table is suitable for usage in all institutions that are directly connected to public power supply, which also supplies buildings built for residential purpose. Quality of supply voltage should be consistent with a typical business or hospital surrounding. Connection to a multi-outlet power strip is specifically prohibited.


## Medical directive

The instrument tables of the ST series are according to the EU regulation 2017/745 (MDR) on medical devices active devices of class I - non-invasive.

Device class according MDR:	I
Basic UDI-DI.:	PP10431TABLE18
Electromagnetic compatibility EMC:	page 59...63

## Transport regulations

- During transportation, please move tables to the lowest position
- Movement of the table without any installed devices
- Unlock all wheels
- After transportation please check mounting of pillar systems and wheels

## Warning and information signs

### Device identification – nameplate

### WARNING

If you do not consider the information marked by this warning, this may result in moderate injury and / or property damage or changes in the product.

### LINE VOLTAGE

If you do not consider the information in a LINE VOLTAGE, this may cause an electrical hazard.

### MAIN VOLTAGE

If you open the housing please separate the unit from the mains voltage.

### CE- marking

This instrument table complies with the requirements of the directive 2017/ 745 medical device regulation (MDR)

Device class according to MDR: I

Basic UDI-DI:

PP10431TABLE18


If you change the product without authorization of the manufacturer, this declaration becomes invalid.


### WEEE- marking

Information on professional waste disposal


### General construction

- 1 tabletop, application part (material HPL laminated plastics)
- 2 handle bar, application part (optional) only for sitting of treatment position (material PA 6 plastics)
- 3 cable guide
- 4 control panel UP/ DOWN
- 5 casing of power supply (optional)
- 6 pillar master/ slave
- 7 power supply with fuse 250VAC, 2x T 10 A/ H  
power cord 2,5m x 1mm<sup>2</sup> and 3pole- IEC connection
- 8 table feet, wheels with fixation

- 9 pillar with network connection CAT.5 E
- 10 4x IEC power outlet, fuse 250VAC, 2x T 6,3 A/ H


power cord, power input ( 7 )  
**In case of error, disconnect here!**

**Please note**

## Operation

The use of the instrument table is reserved for instructed personnel only. Pay attention to the stability of the instrument table.  
Connect the power supply with the provided power cable. The electromotive lifting columns guarantee the height adjustment in a very large area. Adjust the height with the examination equipment installed.

Working time 10%  
1 min ON / 9 min OFF

Operating button - optional (by isolation transformer: page 54)

0 - OFF  
I - ON

## Caution!

Make sure that there are no objects or body parts in the lifting area of the tabletop.  
Note the possible danger of pinching and crushing in the patient environment!  
Please note as an operator that you do not touch the application parts and the patient at the same time!

**Nameplate**

1. Type of device
2. UDI-DI (PPN)
3. Serial number
4. Date of manufacturing
5. Connected load
6. Input mains
7. Mains fuse
8. Weight/ max. load
9. Working time
10. Output, power connection
11. Manufacturer
12. Type

13. Conformity -  
Medical device class I

14. Type of device B

15. Disposal

**Intended use**

The height-adjustable instrument tables serve to inclusion of ophthalmic devices and measuring and testing instruments. Any use other than the one pointed out above is not admissible.

**Devices and accessories**

Ophthalmic examination devices up to max. 80kg (centric load)

## Instructions for assembly and usage

- Warning: To avoid the risk of electric shock, this device may only be connected to a supply network with a protective conductor.
- Follow the accident prevention measures due to legal requirements.
- Compare the data on the type plate and the line voltage.
- The power plug may only be connected to a grounding-type plug.
- Never pull on the mains lead.
- No connection to extension cable or multiple socket
- The instrument tables must assemble in such a way that they can easily be disconnected from the mains.
- The instrument tables may never be assembled and run in rooms with high air humidity. Avoid dripping and splashing water.
- The instrument tables may not be operated in explosive environments.
- Instrument tables may only be assembled and put into operation by qualified personnel who are familiar with the assembly, commissioning and operation of the product. Changes and repairs of the instrument tables may only be performed by our service personnel or authorized distributors. For possible damage caused by non-authorized repairs at the tables, the producer will not be liable. In this case, all warranty claims will expire.
- Please check the possibly necessary regulations and qualification certificates of your country.
- Before first commissioning ensure that the stability of the table is guaranteed.
- Before putting the unit into operation, check the attachment points of the lifting units
- In case of system problems disconnect power cord (7).
- Please observe the instructions for use and instructions of the device manufacturer during installation and commissioning
- You will find necessary spare parts at page 64

## Notes on installation site and accident prevention

The instrument tables may only be moved and set up on a level, hard-hitting ground. The positioning and fixation at the place of destination is carried out by the floor rollers with total lock. These must be adjusted and secured as shown in the following figure to ensure optimum stability of the product.

Risks related to the positioning and the nature of the site can be found in the Risk Management Report and the product FMEA.


When used as intended, the instrument table has a high level of stability. The total locks must be fixed to the floor rollers.

When moving in the transport position, a slope of 10 ° must not be exceeded.

## CAUTION!

Handle bars are used finally finding the treatment position

At high force (for example by pulling or holding) in the direction of the patient is **Risk of Tilt** of the instrument table


**Table wheels** in direction of movement, use fixation of wheels all directions ! (see picture)

## Recommendation

Only use CE-certified products as examination chair or stool, fixable wheels are recommended!

Electrical design

### Device connection IEC socket

Additional multiple sockets or extension cords must not be connected to the instrument table.


The multiple sockets on the instrument table may only be used for the intended ophthalmological devices and a load of the multiple sockets of max. 6.3 A is ensured.

### Change of main fuses

The main fuses are included in the power input socket at the Master pillar. A change may only be made in tension-free operation and by instructed or technically qualified personnel. Please only use the specified fuse sizes!

#### 1. Power input (7)

Fuse holder with clamping plate  
250 VAC Fuse 2x T 10 A/ H

#### 2. Fuse (10)


Turn the fuse holder to Open to the left and reach the device fuse 250VAC 2x T 6,3 A/ H – After this, insert fuse change again and close to the right


### Variations of system tables

The notes on the MED system of Carl Zeiss Meditec AG must be observed!

order-No.	description	dimensions/ weight
ST-100	ST-100, instrument table (e.g. ZEISS Visulas YAG)	1000 x 440 mm Instrument table with handle bars  Install pillar dimensions: 500mm  weight: 35 kg
		
ST-100-A	ST-100-A, instrument table (e.g. ZEISS Atlas 9000)	1000 x 440 mm Instrument table with handle bars, keyboard support  Install pillar dimensions: 500mm  weight: 36,5 kg
		
ST-110-IOL	ST-110-IOL, instrument table (e.g. ZEISS IOL-Master)	1180 x 440mm Instrument table with handle bars, keyboard support  Install pillar dimensions: 500mm  weight: 38 kg
		


order-No.	description	dimensions/ weight
ST-110-VC	ST-110-VC, instrument table (e.g. ZEISS Visucam 200 + 500)	1180 x 440 mm Instrument table with handle bars  Install pillar dimensions: 500mm  weight: 37 kg
		
ST-110-VCs	ST-110-VCs, instrument table (e.g. ZEISS Visucam 224 + 524)	1180 x 440 mm Instrument table, casing for power supply  Install pillar dimensions: 500mm  weight: 40 kg
		
ST-110-KL	ST-110-KL, instrument table (e.g. Visulas 532s + combilaser YAG III/ Visulas 532s)	1180 x 440 mm Instrument table with handle bars, casing for power supply  Install pillar dimensions: 500mm  weight: 46 kg
		

order-No.	description	dimensions/ weight
ST-110-OCT	ST-110-OCT, instrument table (e.g. ZEISS Photo- OCT)	1180 x 440 mm Instrument table with handle bars, casing for power supply  Install pillar dimensions: 500mm  weight: 46 kg
		
ST-115-PR	ST-115-PR, instrument table (e.g. ZEISS PRIMUS OCT)	1280 x 440 mm Instrument table with handle bars, casing for power supply  Install pillar dimensions: 500mm  weight: 48 kg
		
ST-150-OCT	ST-150-OCT, instrument table (e.g. ZEISS: OCT Cirrus)	1180 x 550mm  Install pillar dimensions: 500mm  weight: 38 kg
		

order-No.	description	dimensions/ weight
ST-150T-OCT	ST-150T-OCT, instrument table (e.g. ZEISS OCT Cirrus with isolation transformer 600VA)	1180 x 550 mm Install pillar dimensions: 500mm  weight: 49 kg
		
ST-500-HFA	ST-500-HFA, instrument table (e.g. ZEISS HFA- Perimeter)	1000 x 500 mm Larger pillar (600 mm)  Install pillar dimensions: 600mm  weight: 39 kg
		
ST-500-MP	ST-500-MP, instrument table (e.g. ZEISS Humphrey FDT)	1000 x 500 mm Install pillar dimensions: 500mm  weight: 38 kg
		

order-No.	description	dimensions/ weight
ST-500-MX	ST-500-MX, instrument table (e.g. ZEISSHumphrey Matrix 800)	1000 x 500 mm Keyboard support  Install pillar dimensions: 500mm  weight: 39,5 kg
		
ST-1060	ST-1060, instrument table (e.g. I-Profilier und NCT)	1380 x 600mm  Install pillar dimensions: 500mm  weight: 50 kg
		
ST-2000	ST-2000, instrument table (e.g. 3-4 device)	1700 x 730 mm  Install pillar dimensions: 600mm  weight: 70 kg
		

### Patient surrounding

The patient environment is the area marked above. In this area, the maximum possible protection for the patient is ensured. Dangers can occur in the area of the patient environment. Here is a special duty of care of the attending physician or optician just in connection with the installed devices. In this area, a patient may intentionally or unintentionally have contact with installed components.

- See also: combinations medical and non- medical device (page 53)


## System combination with non-medical devices

To ensure the electrical safety of medical devices and non-medical devices (such as printers), a safe separation according to IEC 60601-1 with 2x MOPP is required. This may require the following measures.

### 1. Isolation transformer (accessory)

To ensure the electrical safety of medical devices and non-medical devices (e.g., printers), operation in medical protective areas is provided with an isolation transformer (accessory).

The isolation transformer used is supplied with mains voltage from the instrument table, complies with the IEC 60601-1 standard and has protection class I.


**Power input:** 1x 230 VAC  
**Power output:** 4x 230 VAC (max. 2,5 A)


### 2. Network isolator (accessory)

The network isolators prevent galvanic connection between the MED and non-MED devices in the signal line (e.g., network feedthrough) as well as the transmission of unwanted voltages and currents. The network isolator is connected directly to the network connection and must comply with the requirements of IEC 60601-1!

**Instrument table with isolation transformer**

- 1 Isolation transformer – 600 VA
- 2 4x IEC outlet 230 VAC – max. 2,5 A
- 3 Power input instrument table (page 39)

- 4 2x fuse holder, 2x T 3,15 A / H
- 5 230 VAC power input isolation transformer
- 6 Operation button ON/ OFF

## Maintenance of device

The instrument tables of the ST series are maintenance-free and designed for a service life of 10 years. Due to their design the instrument tables are mobile! In order to ensure safe and proper operation and a long service life at all times, regular preventive maintenance must be carried out.

The electrical safety of the unit may deteriorate due to aging and use.

At least the following safety checks should be carried out on the system by the manufacturer or qualified persons.

- Check for the presence of the operating instructions
- Visual inspection of device and accessories for damage
- legibility of the inscription
- Leakage current test
- PE conductor test
- Functional and wear check of the castors and brakes
- Functional test of the control element
- Documentation of the results
- Failure to observe the inspection intervals may result in a hazard.

## cleaning of device surface

All table tops have a plastic-coated surface according to EN 14322, are biocompatible according to ISO 10993-1 and are correspondingly resistant. Only the outer surfaces of the instrument tables may be cleaned. A non-dripping cloth with a mild cleaning solution (e.g. soap) is recommended for cleaning. For reasons of hygiene, the instrument table surface and the handles (optional) should be cleaned with alcohol and rapid disinfection after each use!

## Attention!

The use of solvents (ether, acetone), acids or abrasive cleaning agents may cause discoloration or damage.

## Disposal

This symbol is only valid for EU member states.

To avoid any potential negative impact on the environment and possibly on human health, this equipment must be disposed of at a collection point.

Observe the local regulations for the disposal and recycling of hazardous waste in your country.

WEEE-Reg.-No. DE 67707987

## Eliminate electromagnetic interference

The user can reduce electromagnetic interference by observing the recommended minimum distances between portable and mobile RF communications equipment (transmitters) and the equipment. The distance to be maintained depends on the output power of the respective telecommunications equipment, see the following table.

HF-source Wireless Communications facilities	Transm.- frequency (MHz)	Test- frequency (MHz)	max. power P (W)	distance d (m)
Different radio services (TETRA 400)	380-390	385	1,8	0,3
Walkie-talkie (FRS) Rescue service, police, fire department Maintenance (GMRS)	430-470	450	2	0,3
LTE net 13/17	704-787	710/745/ 780	0,2	0,3
GSM 800/900 LTE net 5 Mobile phone CT1+, CT2,CT3	800-960	810/870/ 930	2	0,3
GSM 1800/1900 DECT (mobile phone) LTE Band 1/3/4/25 UMTS	1700-1990	1720/1845/ 1970	2	0,3
Bluetooth, WLAN 802.11b/g/n LTE net 7 RFID 2450 (active & passive transponder & readers)	2044-2570	2450	2	0,3
WLAN 802.11a/n	5100-5800	5240/5500/ 5785	0,2	0,3

The user can take the following measures to eliminate electromagnetic interference

- Increase distance to interference source
- Turn the instrument table or change the angle of the radiation
- Using the instrument table with a different mains connection
- Only use original accessories (e.g. power cord)
- Perform equipotential bonding

**Behavior in case of failure**

If an error occurs, disable the system and use the following troubleshooting list to correct the error. If this does not solve the problem, mark the device as malfunctioning and notify the distributor of your entire system.

**For your  
safety**

**Technical difficulties on the system**

Technical difficulties	Possible cause	Assistance	Reference
No working	Power cord disconnect	Check power cord of incoming voltage	-
	Power failure	Call an electrician	-
	Mains input fuse	Change fuse of mains input	page 45
Devices on Table top - no working	Devices on table top not switched on	Connect device in power outlet or switch on power supply	User manual med- devices
	Fuse of device defective	Check fuse and change	User manual med- devices
	Connection between device and power outlet instrument table failed	Check connection between device and power outlet instrument table	page 45
	Fuse of device defective	Check fuse and replace	page 45
Pillar system – no working	Check – no working		
	Operation button ON/ OFF defective	Check connection between control panel and master pillar	page 40
	Synchron cable defective	Check connection between Master and slave pillar	page 44
Pillar system – no synchron working	Pillar system defective	Contact your supplier	-

## Environmental conditions

Applicable for all variants listed in the instructions for use (page 47-51)

### 1. Put into operation

Attribute	Permissible value range
Temperature	+5°C ... +40°C
Relative humidity	30% ... 75% no condensation
Altitude	2000m above sea level

### 2. Transportation and storage

Attribute	Permissible value range
Temperature	-20°C ... +70°C
Relative humidity	10% ... 90% no condensation
Air pressure	500hPa ... 1060hPa

## Technical data

Attribute	Permissible value range
Nominal voltage	230V AC $\pm$ 10%
Frequency	50 Hz
Medical class	I
Type of device	B
Input pillar	2 A
Power output	6,3 A
Operation duty 10%	1 min ON/ 9 min OFF
Minimum working high	650 mm (pillar: 500mm)
Maximum working high	900 mm (pillar: 500mm)
Lift way	250 mm
Lifting speed	10 mm/s
Weight	35kg – 49 kg (page 47-51)
Lift	max. 80 kg (centric load)

Wagner & Guder Medical GmbH  
Hermstedter Straße 57  
99518 Bad Sulza - Hermstedt, Germany

**Manufacturer's Declaration - Electromagnetic Emissions (Tab. 201 according to DIN EN 60601-1-2)**

The system table is intended for use in an electromagnetic environment as described below. The customer or user of the instrument table should ensure that the device is used in such an environment.

Warning! The use of longer cable length may cause an increased emission or a reduced interference resistance. The use of other sensors or cables except the ones mentioned above is not allowed.

- Cable length incoming voltage (2,5m)

Manufacturer's Declaration - Electromagnetic Emissions		
The system table ST-series is intended for use in an electromagnetic environment as described below. The customer or user of the instrument table should ensure that the device is used in such an environment.		
Emission Measurements	Accordance	Electromagnetic Environment - Guidelines
HF emissions acc. to CISPR11	Group 1	The device uses HF energy exclusively for its internal function. Thus the HF emission is very low and it is unlikely that nearby electronic devices would be disturbed.
HF emissions acc. to CISPR11	Class B	The device is intended for use in all facilities including living quarters and such ones which are connected directly to a public power supply that supplies also buildings used for living purposes.
Emission of overtones acc. to IEC61000-3-2	Class A	
Emission of voltage fluctuation/flicker acc. to IEC61000-3-3	agreed	

**Recommended Safety Distances between portable and mobile HF  
Telecommunication Devices and the refraction unit  
(Tab. 206 according to DIN EN 60601-1-2)**

The system table ST-series is intended for use in an electromagnetic environment with controlled HF disturbances. The user of the device can help to avoid electromagnetic disturbances by keeping the minimum distance between portable and mobile telecommunication devices (transmitters) and the device - depending on the output power of the telecommunication devices as described below.

Nominal power Of transmitter  W	Safety Distance Depending on the Frequency in m		
	150KHz - 80MHz $d=1,2\sqrt{P}$	80MHz - 800MHz $d=1,2\sqrt{P}$	800MHz - 2,7GHz $d=2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum nominal power not mentioned above: To detect the recommended safety distance use the equation in the corresponding column. P is the maximum nominal power of the transmitter in watt (W) according to the specifications of the transmitter manufacturer.

NOTE 1: For 80 Hz and 800 MHz the higher frequency range is valid.

NOTE 2: These guidelines may not be applicable for all cases. The propagation of electromagnetic values is influenced by absorptions and reflections of buildings, objects and people.

<b>Manufacturer's Declaration - Electromagnetic Interference Resistance (Tab. 204 according to DIN EN 60601-1-2)</b>			
The system table ST-series is intended for use in an electromagnetic environment as described below. The customer or user of the table should ensure that the device is used in such an environment.			
<b>Interference Resistance Test</b>	<b>IEC 60601- Testing Level</b>	<b>Accordance Level</b>	<b>Electromagnetic Environment - Guidelines</b>
Conducted HF-disturbances Acc. To IEC 61000-4-6	3Vrms 150kHz bis 80MHz	3Vrms	Portable and mobile radio sets should be used in a no less distances to the device including the cables than it is recommended by the equation for the frequency.  Recommended safety distance: $d=0,35\sqrt{P}$
	6Vrms 150kHz bis 80MHz	6Vrms	
Radiated HF-disturbances Acc. To IEC 61000-4-3	3V/m 80MHz bis 2,5GHz	3V/m	$d=0,35\sqrt{P}$ 80MHz – 800MHz $d=0,35\sqrt{P}$ 800MHz – 2,7GHz P is the nominal power of the transmitter in watt (W) according to the specifications of the transmitter manufacturer; d is the recommended safety distance in meters (m).(a)  The field strength of stationary transmitters should be lower than the accordance level for all frequencies according to a testing on location.(b) Disturbances are possible near devices with the following symbol:
	10V/m 80MHz bis 2,7GHz	10V/m	

NOTE 1: For 80 Hz and 800 MHz the higher frequency range is valid.  
NOTE 2: These guidelines may not be applicable for all cases. The propagation of electromagnetic values is influenced by absorptions and reflections of buildings, objects and people.

<sup>a</sup> The field strength of stationary transmitters such as fixed parts of cellular phones and mobile radio sets, amateur radio stations, AM and FM radio and television cannot be determined exactly in theory. To detect the electromagnetic environment in regard to stationary transmitters a study of the location should be considered. If the measured field strength at the location where the device is being used exceeds the accordance level above the device should be watched to verify the proper functions. If unusual features are watched additional actions might be necessary such as a modified orientation or another location of the device.

<sup>b</sup> For the frequency range of 150 kHz to 80 MHz the field strength should be lower than 3 V/m.

**Manufacturer's Declaration - Electromagnetic Emissions  
(Tab. 202 according to DIN EN 60601-1-2)**

The system table ST-series is intended for use in an electromagnetic environment as described below. The customer or user of the table should ensure that the device is used in such an environment.

Interference Resistance Test	IEC 60601- Testing level	Accordance Level	Electromagnetic Environment - Guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be of wood, concrete or ceramic tiles. If the floor is tiled with synthetic material the relative air humidity must have 30 % at least.
Fast transient electric disturbances / bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	± 2 kV for power lines ± 1 kV for input and output lines	The quality of the supply voltage should conform to a typical business or clinic environment.
Surge voltage acc. to IEC 6100-4-5	± 1 kV normal mode voltage ± 2 kV common mode voltage	± 1 kV normal mode voltage ± 2 kV common mode voltage	The quality of the supply voltage should conform to a typical business or clinic environment.

Interference Resistance Test	IEC 60601- Testing level	Accordance Level	Electromagnetic Environment - Guidelines
Voltage drops, short interruptions and variations in supply voltage acc. to IEC 61000-4-11	<p>&lt; 5% U<sub>T</sub> für 1/2 period (&gt; 95% Break) 40% U<sub>T</sub> für 5 Period (60% Break) 70% U<sub>T</sub> für 25 Period (30% Break) &lt; 5% U<sub>T</sub> für 5 s (&gt; 95% Break)</p>	<p>&lt; 5% U<sub>T</sub> für 1/2 Periode (&gt; 95% Break) 40% U<sub>T</sub> für 5 Period (60% Break) 70% U<sub>T</sub> für 25 Period (30% Break) &lt; 5% U<sub>T</sub> für 5 s (&gt; 95% Break)</p>	<p>Floors should be of wood, concrete or ceramic tiles. If the floor is tiled with synthetic material the relative air humidity must have 30 % at least. The quality of the supply voltage should conform to a typical business or clinic environment.</p>
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields at the supply frequency should conform to the typical values as they occur in the business or clinic environment.
<p>NOTE: U<sub>T</sub> is the AC mains voltage before the use of testing levels</p>			

**Order data - spare parts**

Description	Order- No.
Power cord 2,5m	MC-B1D12500
Device power cord 1m	MG-D1S11000
Fuse T 6,3 A / H	42G1788
Fuse T 10 A / H	42G1794
Control panel UP/ DOWN	DPA-K-06
Handle bars (2x)	GN798
Direction for use system tables	WG-SGA